

MEDICAL IMAGE PROCESSING (67705):


EXERCISE 1-PART 2 - LIVER SEGMENTATION


BY: AMIT HALBREICH, ID: 208917393

<u>CT Case 1 – Liver Segmentation Result</u>









Ground Truth File:

MIP Data/Case1_liver_segmentation.nii.gz

Calculated Liver Segmentation:


MIP Data/Segmentations/Case1_liver_seg.nii.gz


DICE Score & VOD Score


DICE = 0.9101813896660667


VOD = 0.16483221975152706


<u>CT Case 2 – Liver Segmentation Result</u>















<u>CT Case 3 — Liver Segmentation Result</u>













<u>CT Case 4 – Liver Segmentation Result</u>













<u>CT Case 5 – Liver Segmentation Result</u>













<u>CT Hard Case 1 – Liver Segmentation Result</u>









Ground Truth File:

MIP Data/HardCase1_liver_segmentation.nii.gz

Calculated Liver Segmentation:

 ${\bf MIP\ Data/Segmentations/HardCase1_post_process_liver_segmentation.nii.gz}$


DICE Score & VOD Score


DICE = 0.8804743206954235


VOD = 0.21352914276843227


<u>CT Hard Case 2 - Liver Segmentation Result</u>

As we can spot the cleaning process was a little harder for this case – even after taking "tight" boundaries for the liver ROI. Some of the skin and other organs were added maybe because they are somewhat close in their HU Unit value.













<u>CT Hard Case 3 — Liver Segmentation Result</u>














<u>CT Hard Case 4 – Liver Segmentation Result</u>















<u>CT Hard Case 5 – Liver Segmentation Result</u>

As we can spot because there are only 22 axial slices in this CT scan it has very low resolution but the shape of the Liver Segmentation resembles the triangular anatomical shape of the liver. And the axial slice of the segmentation on the top-right image resembles an actual liver axial slice.









Functions & Explanations

Code Style & Insights

Use of Vectorial and Matrixial Coding:

<u>NumPy Arrays</u>: The code utilizes NumPy arrays for efficient handling of multidimensional data, enabling vectorized operations.

<u>Vectorized Operations</u>: Vectorized operations are extensively used throughout the code for tasks like thresholding, finding voxel neighbors, computing gradient magnitudes, running Multi Seeded Region Growing algorithm with matrixial and vectorized data – taking advantage of taking all seeds and their neighbors to the queue and updating it in each iteration by taking the new front voxels that hasn't been visited yet – all this for enhancing performance by avoiding explicit loops.

<u>Matrix Operations:</u> Matrices are used to represent voxel coordinates and masks, facilitating efficient operations like element-wise multiplication and logical indexing.

<u>Optimized Algorithms:</u> The MSRG algorithm is implemented efficiently using matrix operations and vectorized logic, reducing redundant computations and improving runtime performance.

Results and Performance:

<u>Segmentation Accuracy:</u> The segmentation accuracy is evaluated using DICE and VOD scores, comparing the segmented liver with ground truth data. According to the Formula:

$$\begin{split} DICE &= \frac{2 \cdot |VOL(GT_Segmentation) \cap VOL(Estimated_Segmentation)|}{|VOL(GT_Segmentation)| + |VOL(Estimated_Segmentation)|} \\ VOD &= 1 - \frac{|VOL(GT_Segmentation) \cap VOL(Estimated_Segmentation)|}{|VOL(GT_{Segmentation}) \cup VOL(Estimated_Segmentation)|} \end{split}$$

<u>Runtimes:</u> The emphasis on vectorial and matrixial coding helps in achieving shorter runtimes, especially for computationally intensive tasks like MSRG, ensuring reasonable performance even for large nifti files.

<u>Scalability:</u> The code design allows for scalability, enabling it to handle different CT scan sizes efficiently.

<u>Output:</u> The code saves segmented liver data and evaluation results, providing insights into segmentation quality and aiding further analysis of different stages in the liver segmentation process.

Imported Library Functions

```
import os
import random

import numpy as np
import nibabel as nib
from skimage.morphology import binary_opening, binary_closing, ball,
label
from skimage.measure import label
```

- I used skimage.morphology .binary_opening and skimage.morphology.binary_closing for operations for post_processing or preprocessing data before performing binary threshold or get the largest Connected Component. skimage.morphology.ball is used as binary opening and binary closing object to dilate or erose.
- I used numpy library for operations and calculations.
- I used nibabel library in order to load, save and read nifty files.
- I used skimage.measure.label library function in order to get labeled scans and get the number of Connected Components to make sure I get 1 CC after body isolation process etc.

Nifti files Utils - Helper Functions

```
:return: The image data, the image filename, and the Nifti file
   :param img data: The image data to be saved.
   :param ct filename path: The path to the original CT filename.
def save scan data to path(img data, output path, file name):
   :param img data: The scan data to be saved.
```

Isolate Body

```
Function Name: def isolate_body(ct_data, body_lower_threshold=BODY_LOWER_THRESHOLD, body_upper_thresh=BODY_UPPER_THRESHOLD, ball_size=2, ct_path=''):

"""

This function isolates the body from CT data and separates between the patient's body mask and the bed of the CT chamber. At the end the noise

from the image is cleaned.

:param ct_data: The CT data.

:param body_lower_threshold: The lower threshold for the body.

:param body_upper_thresh: The upper threshold for the body.

:param ball_size: The size of the ball structuring element.


:param min_size_obj_range: The range of objects to remove.

:return: The segmented isolated body.

"""
```

Input: ct_scan (Hard Case 1-5 or Case 1-5)

Output: body_mask with patient's body separated from the bed (for Hard Case 1-5 or Case 1-5) – We can spot the input images has the patients' beds in it and after body isolation we remain only with body segmentation without the bed or noise.



Create Liver ROI


```
Function Name: def create_liver_roi(ct_data, aorta_data, liver_lower_th=LIVER_MIN_THRESHOLD, liver_upper_th=LIVER_MAX_THRESHOLD, liver_roi_x=(-30, 195), liver_roi_y=(-100, 100), liver_roi_z=(-75, 50), ct_path=''):

"""

This function creates a region of interest (ROI) for the liver.

:param ct_data: The CT data.

:param aorta_data: The aorta data.
```



Find Seeds



Multiple Seeds Region Growing

```
Function Name: def multiple_seeds_RG(ct_data, roi_data,
max_iter=MAX_ITERS, ct_path=''):
    """
    Perform Multi-Seeded Region Growing algorithm fully with all
parts.


Args:
    ct_data (ndarray): The CT data.
    roi_data (ndarray): The ROI data.
    max_iter (int): Maximum number of iterations.
    ct_path (str): The original CT Scan path.

Returns:
    ndarray: Segmented liver.
"""
```

Input: ct_data, roi_data, max_iter(Optional Variable) (Hard Case 1-5 or Case 1-5)

Output: seeds Liver Segmentation estimations after running MSRG algorithm – uses vectorial\matrixial code in order to run efficiently and dilate the liver segmentation mask according to BFS Algorithm that adds voxels inside the ROI with vectorized code that an homogeneity condition holds for them: Their HU Value is close according to Homogeneity Function as follows:

$$\delta(vox) = \frac{HU_Value(vox) - MEAN_{y \in R_{seg_{mask}}} \{HU_Value(y)\}}{STD_{y \in R_{seg_{mask}}} \{HU_Value(y)\}}$$



<u>Evaluate Segmentation - Calculate DICE & VOD Scores</u>

```
Function Name: def evaluate segmentation (ground truth seg,
processed liver seg):
     :param ground_truth_seg: The ground truth segmentation.
:param processed_liver_seg: The processed liver segmentation.
:return: The DICE score and VOD.
            Input: ground_truth_seg, processed_liver_seg (HardCase 1 or Case 1)
 Output: DICE Score, VOD of the comparison between Ground Truth Liver Segmentation &
          Estimated Calculated Liver Segmentation from MSRG algorithm results.
         Input Images:
                                   Outputs
                                                         Input Images:
                                                                                     Outputs
       Ground Truth File:
                                                      Ground Truth File:
Case1 liver segmentation.nii.gz
                                                          HardCase1
                                                   liver segmentation.nii.gz
                                   DICE
                                                                                     DICE
                                   = 0.91
                                                                                     = 0.88
   Estimated Segmentation:
                                                   Estimated Segmentation:
                                   VOD
                                                                                     VOD
      Case1_post_process
                                   = 0.165
                                                   HardCase1_post_process
                                                                                     = 0.21
   _liver_segmentation.nii.gz
                                                   _liver_segmentation.nii.gz
```

Multiple Seeds Region Growing

```
Function Name: def segment liver(ct filename, aorta seg filename,
   Input: ct_filename, aorta_seg_filename, output_filename (Hard Case 1-5 or Case 1-5)
   Output: Returns None but runs full Liver Segmentation with all stages in the process:
                                     1. isolate body
                                   2. create liver roi
                            3. find seed + multiple seeds RG
                    evaluate segmentation – Calculate DICE & VOD Scores.
                      Output Image:
  Input Image:
                                             Input Image:
                                                                    Output Image:
    CT Case 1
                         CT Case 1
                                               CT Case 2
                                                                      CT Case 2
                    Liver Segmentation
                                                                 Liver Segmentation
```

Helper Functions – for isolate body and create liver roi functions:

```
def get_segmentation_bbox(segmentation):
    """
    Retrieve the bounding box of a segmentation.
    Parameters:
        segmentation (ndarray): The segmentation.
    Returns:
        tuple: The bounding box coordinates in the order
        (x_min, x_max, y_min, y_max, z_min, z_max).
    """

# Find the minimum and maximum indices for each axis
    x_idx = np.any(segmentation, axis=(1, 2)).nonzero()[0]
    y_idx = np.any(segmentation, axis=(0, 2)).nonzero()[0]
    z_idx = np.any(segmentation, axis=(0, 1)).nonzero()[0]
    x_min, x_max = x_idx[0], x_idx[-1]
    y_min, y_max = y_idx[0], y_idx[-1]
    z_min, z_max = z_idx[0], z_idx[-1]
```

```
# Return the values as six separate variables
return x min, x max, y min, y max, z min, z max
```

Stage 1 : Isolate Body

```
get largest component(mask scan data):
    largest component mask = labeled scan == largest label
    return largest component mask
def isolate body(ct data, body lower threshold=BODY LOWER THRESHOLD,
    :param ct data: The CT data.
    :return: The segmented isolated body.
body upper thresh)
                                                  f'{SEGMENT FOLDER}/'
                                                      ball size)
def post process bin ct(binary img, ball size=10):
    :param binary_img: The binary CT image data.
:param ball size: The size of the ball structure.
```

Stage 2: Create Liver ROI

```
def create liver roi(ct data, aorta data,
liver upper th=LIVER MAX_THRESHOLD,
                                            f'{SEGMENT FOLDER}/'
perform binary thresholding (body without zeros,
                                             f'{SEGMENT FOLDER}/'
```

```
aorta data.shape[0], x max + liver roi x[1])
```

$Stage\ 3A:\ Find\ Suitable\ in-region\ Seeds$

```
seeds_in_roi = np.zeros like(liver roi).astype(bool)
def get seeds lowest gradients sorted(roi data, seeds):
   :param roi_data: The region of interest data.
   :param seeds: The seed points.
gradients magnitudes.shape)]
def sample random points(roi data, num seeds):
```

```
This function samples random points from the ROI.

:param roi_data: The region of interest data.

:param num_seeds: The number of seeds to sample.

:return: The sampled points.

"""

# Find the indices of the ROI points

roi_points_indices = np.transpose(np.where(roi_data == 1))

# Sample random indices without replacement

random_indices = random.sample(range(len(roi_points_indices)),

num_seeds)

# Extract the sampled points from the ROI points array

sampled_points = [tuple(roi_points_indices[idx]) for idx in

random_indices]

return sampled_points

def find_seeds(roi_data, num_seeds=NUM_SEEDS):

"""

This function finds seeds for region growing.

:param roi_data: The region of interest data.

:param num_seeds: The number of seeds to find.

:return: The seeds.

"""

# Sample random points from the ROI
seeds_list = sample_random_points(roi_data, num_seeds)
return seeds_list
```

Stage 3B: run MSRG Algorithm

```
def get_all_next_front_voxels(voxels_coords):
    """
    Get all neighboring voxels of given coordinates.

Args:
    voxels_coords (ndarray): The coordinates of the voxels.

Returns:
    ndarray: All neighboring voxels.

"""

# Define the neighbors_idx_mask for 26 all_next_front_voxels neighbors_idx_mask = np.array([
        [-1, -1, -1], [-1, -1, 0], [-1, -1, 1],
        [-1, 0, -1], [-1, 0, 0], [-1, 0, 1],
        [0, -1, -1], [0, -1, 0], [0, -1, 1],
        [0, 0, -1], [0, 0, 1],
        [0, 0, -1], [0, 0, 1],
        [0, 1, -1], [0, 1, 0], [0, 1, 1],
        [1, 0, -1], [1, 0, 0], [1, 0, 1],
        [1, 0, -1], [1, 1, 0], [1, 1, 1]
])

# Add neighbors_idx_mask to each point
all_next_front_voxels = voxels_coords[:, None, :] +
neighbors_idx_mask
# Exclude the original points
return np.unique(all_next_front_voxels.reshape(-1, 3), axis=0)
```

```
:param voxel coords: The coordinates of the voxels.
   :return: The neighbors of the voxels.
1])).T.reshape(-1, 3)
def homogeneity condition(ct data, segmented liver, voxels to check,
```

```
seeds = np.array(find_seeds(seeds roi))
   save nifti(liver segmentation.astype(np.uint8), ct path,
                                          f'{SEGMENT FOLDER}/'
       MIN QUEUE SIZE = 0
   while iter counter < max iter and queue.shape[0] >
liver segmentation
```

```
def get valid coords and homogeneity mask(ct data,
~visited matrix[tuple(valid voxels.T)]
```

```
def post process liver segmentation(liver segmentation):
   :param x_min: The minimum X-axis value of the ROI.
   :param x max: The maximum X-axis value of the ROI.
   :param y min: The minimum Y-axis value of the ROI.
   :param y max: The maximum Y-axis value of the ROI.
   :param z min: The minimum Z-axis value of the ROI.
   :param z max: The maximum Z-axis value of the ROI.
```

Stage 4: Evaluate Segmentation VOD & DICE Scores

```
union_vol = np.sum(np.logical_or(ground_truth_seg,
processed_liver_seg))

# Compute DICE & VOD Coefficients
   dice = 2.0 * intersection_vol / (ground_truth_vol +
processed_liver_vol)
   vod = 1 - intersection_vol / union_vol

return dice, vod
```

<u>Stage 5 : Segment Liver - Run Full Segmentation Algorithm</u>

```
segment liver(ct filename, aorta seg filename, output filename):
                                             f'/{ct filename}')
```

```
full_output_path = liver_seg_output_path
    with open(full_output_path.replace('.nii.gz', '.txt'), 'w+')
as f:
    f.write(f" Ground Truth File: {liver_ground_truth_path}\n
VS "

    f"\nCalculated "
    f"Liver Segmentation: {liver_seg_output_path}\n"
    f" DICE Score & VOD Score\n")
    f.write(f"DICE is: {dice}\nVOD is: {vod}\n")
```